

CAREVEAM

308 EXCIMER UVB

INNOVATIVE UVB PHOTOTHERAPY SYSTEM FOR SKIN DISEASE

Lameditech's innovative phototherapy system "CAREVEAM" is 308 excimer system for the treatment of various autoimmune and inflammatory skin disease, which is small, handy treatment solution intense monochromatic UVB phototherapy



EX-3000


LAMEDITECH
www.lameditech.com

Innovative UVB phototherapy system

- Highly effective for various type of autoimmune and inflammatory dermatitis
- Safe and 100% pain-free treatment
- User friendly design and easy to operate
- Low acquisition and operating costs, maintenance free technology(no consumables)

Vitiligo


Eczema

Psoriasis

Atopic Dermatitis

Alopecia Areata

Vitiligo case



Atopic dermatitis case



Product Specification

• Light source	XeCl Excimer Light	• Display	4 Digit segment
• Wavelength	308 nm	• AC power	220~230 V / 1A
• Level	Level 1 ~ Level 5	• Device weight	4 kg
• Power density	10 ~ 30 mW/cm ²	• Device size	300 x 250 x 210 mm
• Single irradiation time	1 ~ 180 Sec	• Spot size	21 x 45 mm



LAMEDITECH CO., Ltd.

#1002, ACE High-End 6th,
Beotkkot-ro 234, Geumcheon-gu,
Seoul, South Korea 08513
Tel. +82-2-852-7982
Fax. +82-2-852-7983
E-mail. sales@lameditech.com

www.lameditech.com

LAMEDITECH



Ministry of Food and
Drug Safety



LAMEDITECH has been developing and supplying the specialized laser devices for the healthy and beautiful life of human and animal based on laser miniaturization technology. Through the continuous research and development, LAMEDITECH commits to become a global top laser expertise company.

EC REP OBELIS S.A.

Bd. Général Wahis, 53-1030
Brussels, Belgium
Tel. +32.2.732.59.54
Fax. +32.2.732.60.03
E-mail. mail@obelis.net

308 EXCIMER SYSTEM
CARREVEAM